



eXome Whole  
Sequencing

# WHOLE EXOME SEQUENCING


Library Preparation Assays

# G2M Whole Exome Sequencing Library Preparation Assay

Exons represent just 1–2% of the genome, they harbour up to 85% of disease-causing mutations, making WES a highly effective and efficient approach for identifying genetic variants associated with disease. It enables high-throughput detection of single nucleotide variants (SNVs), insertions/deletions (in-dels), and other coding-region alterations that underlie a wide range of genetic disorders. For many patients—especially those with rare or complex genetic conditions—the journey to a diagnosis can be long, expensive, and emotionally exhausting.

Conventional exome sequencing panels often require a trade-off between comprehensive genomic coverage and assay performance, leading to uneven coverage, higher duplication rates, ultimately impacting the sensitivity and reliability of variant detection.

Developed with a deep understanding of both clinical and research needs, G2M's Whole exome sequencing assay is designed with exceptional uniformity and high on-target efficiency with the panel content aligned with the latest curated genomic data for enhanced clinical relevance. The panel encompasses ~21,500 genes catering to various hereditary conditions and germline cancers.



## Key Features



### Comprehensive Genomic Content

- Covers >99% of clinically relevant protein-coding regions at >20X
- Includes updated gene content based on RefSeq, GENCODE, OMIM, ClinVar & HGMD



### Exceptional Coverage Uniformity

- Minimizes dropouts in GC-rich or hard-to-capture exons
- Enables confident detection of heterozygous and low-frequency variants



### High On-Target Rate

- Maximized hybridization efficiency ensures a high percentage of reads align to intended targets
- Reduces off-target noise, lowering sequencing waste and cost



### Clinical-Grade Performance

- Validated for SNVs, indels, and CNV calling
- Applications in rare disease, oncology, reproductive health, and inherited disorders



### Flexible and Scalable Workflow

- Compatible with major sequencing platforms such as Illumina, MGI, Element Biosciences (AVITI)


## Key Application Areas




## Performance Specifications

Metric	Performance
On-Target Rate	>85%
Coverage Uniformity ( $\geq 0.2 \times$ mean)	>95% of targets
Mean Coverage at 100× Sequencing	>98% of targets covered at $\geq 20 \times$
Duplication Rate	<10%
Fold 80 Base Penalty	1.25
Mapped Reads	>99%
Sensitivity for SNVs	>99.5%
Sensitivity for InDels	>98%
GC Bias Performance	Optimized for high- and low-GC regions
Panel Update Cycle	Biannually (aligned with RefSeq, ClinVar updates)
Validated platforms	Illumina, MGI, Element Biosciences (AVITI)

## Quality Metrics




**Fig 01: NF1 Gene Coverage Sequencing for Precise Variant Detection** : High-depth sequencing of the NF1 gene demonstrates exceptional coverage, enabling precise detection of both pathogenic variants and variants of uncertain significance (VUS), including clinically relevant deletions. This comprehensive analysis supports accurate diagnosis of Neurofibromatosis Type 1 and related syndromes such as Type 1 & 2 Classic Ehlers-Danlos syndrome.




**Fig 02: Comprehensive Gene Coverage Across Key Databases :** The chart illustrates the high percentage of genes mapped to major clinical and phenotype databases - OMIM, ClinVar, and HPO ensuring robust integration of genetic information. Nearly 100% coverage in OMIM and HPO, along with over 90% in ClinVar, highlights the robustness of these databases for accurate gene interpretation and evidence based clinical insights.

## Post-Sequencing Data Quality



**Fig 03: G2M Delivers Superior GC Coverage Uniformity Across Leading Platforms :** GC content distribution across target regions for Company A (red), G2M (blue), and Company B (green) platforms reveals superior uniformity in G2M's coverage across varying GC percentages. This indicates reduced GC bias and more consistent hybridization efficiency, supporting reliable variant detection across challenging genomic landscapes.






**Fig 04: Comprehensive Gene Coverage Across Key Databases :** The chart illustrates the high percentage of genes mapped to major clinical and phenotype databases - OMIM, ClinVar, and HPO ensuring robust integration of genetic information. Nearly 100% coverage in OMIM and HPO, along with over 90% in ClinVar, highlights the robustness of these databases for accurate gene interpretation and evidence based clinical insights.

## Multi-Metric Comparison of Exome Performance




**Fig 05: Exceptional Coverage Uniformity Validated by Low Fold 80 Metric :** Fold 80 base penalty is the fold-over-mean coverage required to ensure that 80% of targeted bases are covered at or above the mean coverage depth. It is a key metric used to evaluate the uniformity of coverage in next-generation sequencing (NGS). A lower value indicates more uniform coverage—which means less over-sequencing is needed to achieve complete and reliable results. G2M shows <1.25 Fold 80 penalty value implying good capture design and hybridization efficiency than the competition.



**Fig 06: Efficient Whole Exome Sequencing with Ultra-Low Depth Coverage :** G2M WES assay shows >98% coverage at a very low depth (1X)


**Fig 07: High On-Target Alignment Across Different Cancer Patient Samples :** G2M WES assay shows an average on target ratio of 85% for both clinical and control samples (NA12878)



## Panel Specifications

No of genes	~21,500
Inclusions	SNVs, Indels, CNVs, Mitochondrial genome
Target size	38.2 Mb
Probe region	49.8 Mb
Recommended coverage	150-180X
Data required (per sample)	6 GB
Pre-capture pooling	8
Methodology	Hybridization capture based target enrichment
Recommended sequencing mode	Paired-end 150 (PE 150)
Sample types	Blood, Saliva, Amniotic fluid, Chorionic Villus, Tissue
Starting sample input	50-500 ng Genomic DNA
Average library insert size	~300 bp

## End-to-End G2M Workflow




\* CSM = Clinically Significant Mutations

## G2M WES Assay Library Construction Workflow


### The library preparation is based on Hybridization capture based target enrichment

Our Whole Exome Sequencing (WES) panels utilize a robust hybridization capture-based target enrichment approach to ensure high specificity and uniform coverage across the exome. This workflow enables the selective enrichment of coding regions from genomic DNA, allowing for comprehensive detection of single nucleotide variants (SNVs), insertions/deletions (indels), and other clinically relevant mutations. The method combines proven biochemical precision with optimized probe design, providing a reliable foundation for accurate and high-throughput sequencing results.

#### Genomic DNA Lib Preparation




#### Hybridisation and Target Enrichment



## Interpret and report relevant variants with Cliseq Interpreter Platform

The NGS data analysis is supported by combining guideline recommended variants with the analytical capability of G2M's Cliseq Interpreter Platform. Cliseq Interpreter is a cloud based NGS data analysis software which offers an unparalleled platform performance designed to streamline and enhance the interpretation of complex biological data. Once Quality Check, Alignment, Variant calling, and annotations are achieved, the annotated VCF files and clinically significant mutations (CSM) report will be available to download.




\*As per allele frequency & phenotypic indicator.

## CSM Report

**Patient Details**

Accession ID:	538
Sample Type:	Blood
Patient's Name:	NF29-Mohd-TEST
Gender/Age:	Male/



WHOLE EXOME SEQUENCING PANEL

**CLINICAL INDICATIONS / PHENOTYPE:**

Neurofibromatosis

**Indication for sequencing**

Data Uploaded On: 27/12/2023

Report Generated On: 27/12/2023

**RESULTS:**

**List of clinically significant mutations (CSM)**

GENE	VARIANT	ZYGOSITY	DISEASE	INHERITANCE	INTERPRETATION
NF1 (NM_001042492.3)	Exon 15 c.1655T>C (p.Leu552Pro)	Heterozygous	Neurofibromatosis- Noonan syndrome (OMIM#601321) Juvenile myelomonocytic leukemia (OMIM#607785) Neurofibromatosis type 1 (OMIM#162200) Neurofibromatosis, familial spinal (OMIM#162210) Watson syndrome (OMIM#193520)	AD, SMu	Pathogenic

**Clinically Significant Mutations**

- Exon number
- Nucleotide change identified at the specific locus
- Amino Acid substitution due to the SNV

**Interpretation derived from the CSMs**

**Disorders linked to the reported CSMs**

**Details of the identified CSM highlighting –**

- Genomic location of mutation
- Type of mutation
- Mutation depth
- Frequency of variation
- Prediction tools employed in analysis
- Clinical evidence of identified mutation in the phenotype

**INTERPRETATION:**

**NF1;c.1655T>C(p.Leu552Pro)**

Genomic location: chr17:31221863 T>C

Variant type: Missense

Depth: total 45; alt 22 (49%)

Frequency: The variant is very rare, absent from the gnomAD reference population dataset.


Prediction tools: Multiple lines of computational tools predict a deleterious effect on the gene or gene product.

Clinical evidence: This variant has previously been described in ClinVar (VCV996460) with the following classifications: VUS (1) / LP (1).

\*Genetic test results are reported based on the recommendations of American College of Medical Genetics.  
\*\*Sanger validation is recommended for the above listed variants.  
\*\*\*The significance/classification of the variant(s) may change based on genetic testing in the parents and other family members.  
\*\*\*\*Data reevaluation performed upon the up gradation of databases used and results may vary in accordance.

*This report is automatically generated by Cliseq Interpreter*

Page 1 of 5



Scan for WES Sample Report

Automate your  
**NGS Libraries** with

**EZY-AutoPrep - 48**

Construct upto 48 sample libraries in one run



Construct  
upto 48 Sample  
Libraries in  
one run

## Ordering Information

Commercial Name	Cat No.
Clinical Exome Sequencing Expanded Panel (Whole Exome Sequencing)	G2MCES07001(WES)-ill; G2MCES07001(WES)-MG, G2MCES07001(WES)-TF



### Genes2Me Pvt. Ltd.

1105, 11th Floor, Tower B, SAS Tower, Medicity,  
Sector - 38, Gurgaon - 122 001, Haryana, India

Tel : + 91 18001 214030 / + 91 88000 23600  
E-mail : [contact@genes2me.com](mailto:contact@genes2me.com)

[www.genes2me.com](http://www.genes2me.com)

