


For **Extraction**,
Enrichment and
Purification of RNA in
Clinical Specimens

MagRNA-II Viral RNA Extraction Kit

- Sample Type : Human Serum, Nasopharyngeal Swabs, Oropharyngeal Swabs, Bronchoalveolar Lavage Fluids.
- Magnetic beads based technology allowing recovery of >95-100% viral RNA.
- Carrier RNA included to enhance quality of eluted viral RNA.
- Extraction Time of ~ 30 - 35 minutes per batch of 24-48 samples.
- Compatible with Manual and Multiple Automated systems such as MGISP-100B, MGISP-960, KingFisher Flex, ILS GENFast Automated System, Agilent Bravo, Hamilton Automated Liquid Handlers, Tecan Freedom Evo100.
- **Plastic ware for specific automation instruments also available**

Ordering Information	Cat. No. G2M030320	200 Test
	Cat. No. G2M030420	500 Test
	Cat. No. G2M030620	960 Test



GENES2me
ADVANCED MOLECULAR DIAGNOSTIC SOLUTIONS

MANF. LICENSE BY
CDSCO

VALIDATED BY
ICMR

CE **IVD**



ISO
13485:2016

ISO
9001:2015



+91 98990 62009

covid@genes2me.com

www.genes2me.com

Gurgaon, Haryana - 122 001

+ 91 88000 23600

covid@genes2me.com