


For **Extraction,**
Enrichment and
Purification of RNA in
Clinical Specimens



SpinRNA Viral RNA Extraction Kit


- Sample Type : Human Plasma, Serum, Nasopharyngeal Swabs, Oropharyngeal Swabs, Bronchoalveolar Lavage Fluids.
- Silica membrane column based technology allowing recovery of >90-95% viral RNA.
- Carrier RNA included to enhance quality of eluted viral RNA.
- Works well with minimal sample volume of 200 µl and elution of high yield/quality of RNA in 30-50 µl.
- Extraction Time of 35-45 minutes per batch.
- Compatible with both Manual and Automated platforms.



Ordering Information

Cat No.
G2M080820

Pack Size
250



GENES2me
ADVANCED MOLECULAR DIAGNOSTIC SOLUTIONS

MANF. LICENSE BY
CDSCO

VALIDATED BY
ICMR

CE IVD



ISO
13485:2016

ISO
9001:2015



+91 98990 62009

covid@genes2me.com

www.genes2me.com

Gurgaon, Haryana - 122 001

+ 91 88000 23600

covid@genes2me.com