


Report Title

Test Name:	TB-Q Comp	Test Sample:	328042025
Test Time:	2025-04-28-16-09-45_1	Report Time:	2025-04-28

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	katG	Unknown	25.97	NA	katG Positive	katG Positive, RIF-2 Positive, RIF-3 Positive, mpt/IS611 Positive	katG Positive, RIF-2 Positive, RIF-3 Positive, mpt/IS611 Positive, RIF-1 Positive, RIF-4 Positive, INH-inhA Positive, RNaseP Positive
	VIC	RIF-2	Unknown	23.43	NA	RIF-2 Positive		
	ROX	RIF-3	Unknown	23.12	NA	RIF-3 Positive		
	CY5	mpt/IS611	Unknown	25.15	NA	mpt/IS611 Positive		
Well B	FAM	RIF-1	Unknown	26.42	NA	RIF-1 Positive	RIF-1 Positive, RIF-4 Positive, INH-inhA Positive, RNaseP Positive	Positive, katG Positive, RIF-2 Positive, RIF-3 Positive, mpt/IS611 Positive, RIF-1 Positive, RIF-4 Positive, INH-inhA Positive, RNaseP Positive
	VIC	RIF-4	Unknown	24.13	NA	RIF-4 Positive		
	ROX	INH-inhA	Unknown	24.95	NA	INH-inhA Positive		
	CY5	RNaseP	Unknown	19.47	NA	RNaseP Positive		
Well C	FAM	katG	Unknown	26.72	NA	katG Positive	katG Positive, RIF-2 Positive, RIF-3 Positive, mpt/IS611 Positive	Positive, katG Positive, RIF-2 Positive, RIF-3 Positive, mpt/IS611 Positive, RIF-1 Positive, RIF-4 Positive, INH-inhA Positive, RNaseP Positive
	VIC	RIF-2	Unknown	24.3	NA	RIF-2 Positive		
	ROX	RIF-3	Unknown	23.43	NA	RIF-3 Positive		
	CY5	mpt/IS611	Unknown	25.71	NA	mpt/IS611 Positive		
Well D	FAM	RIF-1	Unknown	26.81	NA	RIF-1 Positive	RIF-1 Positive, RIF-4 Positive, INH-inhA Positive, RNaseP Positive	Positive, INH-inhA Positive, RNaseP Positive
	VIC	RIF-4	Unknown	24.71	NA	RIF-4 Positive		
	ROX	INH-inhA	Unknown	25.25	NA	INH-inhA Positive		
	CY5	RNaseP	Unknown	19.63	NA	RNaseP Positive		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.