



ViralCollect Transport Kit

Ready to use Viral transport kit intended for the collection and transport of clinical specimens containing viruses from the collection site to the testing laboratory.

- 100 sterile specimen collection polyester swabs and 50 tubes of Viral Transport Medium.
- 3ml VTM in 15 ml centrifuge tube.
- VTM Contains protective antibiotics, antifungal agents to control microbial contamination and buffer to stabilize the pH.
- pH of 7.3 + 0.3 at 25°C along with Phenol Red as pH indicator.
- Contains cryo-protectant to preserve for prolonged storage.
- Osmolality in mOsm/Kg H₂O: 500 - 600 mOsm.



MANF. LICENSE BY
CDSCO

VALIDATED BY
ICMR

CE IVD



ISO
13485:2016

ISO
9001:2015

+91 98990 62009

covid@genes2me.com

www.genes2me.com

Gurgaon, Haryana - 122 001

91 88000 23600

covid@genes2me.com