


Report Title

Test Name:	AMR-Q	Test Sample:	11082025AMR
Test Time:	2025-08-11-12-47-01_1	Report Time:	2025-08-11

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	VIM3	Unknown	22.31	NA	VIM3 Positive	VIM3 Positive, OXA48 Positive, IMP2 Positive, NDM2 Positive, CY5.5 Positive	VIM3 Positive, OXA48 Positive, IMP2 Positive, NDM2 Positive, CY5.5 Positive
	VIC	OXA48	Unknown	22.45	NA	OXA48 Positive		
	ROX	IMP2	Unknown	24.48	NA	IMP2 Positive		
	CY5	NDM2	Unknown	22.39	NA	NDM2 Positive		
	CY5.5	CY5.5	Unknown	24.46	NA	CY5.5 Positive		
Well B	FAM	Mec-A	Unknown	21.79	NA	Mec-A Positive	Mec-A Positive, CTX-M Positive, Van A Positive, MCR-1 Positive, CY5.5 Positive	Mec-A Positive, CTX-M Positive, Van A Positive, MCR-1 Positive, CY5.5 Positive
	VIC	CTX-M	Unknown	23.33	NA	CTX-M Positive		
	ROX	Van A	Unknown	20.31	NA	Van A Positive		
	CY5	MCR-1	Unknown	22.66	NA	MCR-1 Positive		
	CY5.5	CY5.5	Unknown	22.8	NA	CY5.5 Positive		
Well C	FAM	TEM	Unknown	21.47	NA	TEM Positive	TEM Positive, Van B Positive, SHV Positive, CiTm Positive, CY5.5 Positive	TEM Positive, Van B Positive, SHV Positive, CiTm Positive, CY5.5 Positive
	VIC	Van B	Unknown	21.95	NA	Van B Positive		
	ROX	SHV	Unknown	22.33	NA	SHV Positive		
	CY5	CiTm	Unknown	21.41	NA	CiTm Positive		
	CY5.5	CY5.5	Unknown	21.12	NA	CY5.5 Positive		
Well D	FAM	dfrA1	Unknown	21.37	NA	dfrA1 Positive	dfrA1 Positive, dfrA2 Positive, Sul 2 Positive, PMQR Positive, CY5.5 Positive	dfrA1 Positive, dfrA2 Positive, Sul 2 Positive, PMQR Positive, CY5.5 Positive
	VIC	dfrA2	Unknown	21.43	NA	dfrA2 Positive		
	ROX	Sul 2	Unknown	22.03	NA	Sul 2 Positive		
	CY5	PMQR	Unknown	21.81	NA	PMQR Positive		
	CY5.5	CY5.5	Unknown	21.48	NA	CY5.5 Positive		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.