


Report Title

Test Name:	BCR ABL	Test Sample:	05040120251
Test Time:	2025-01-04-18-16-39_1	Report Time:	2025-01-04

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	BCR-ABL1	Unknown	25.07	NA	BCR-ABL1 Positive	BCR-ABL1 Positive, ABL1 Positive	BCR-ABL1 Positive, ABL1 Positive
	ROX	ABL1	Unknown	26.48	NA	ABL1 Positive		
Well B	FAM	BCR-ABL1	Unknown	26.16	NA	BCR-ABL1 Positive	BCR-ABL1 Positive, ABL1 Positive	BCR-ABL1 Positive, ABL1 Positive
	ROX	ABL1	Unknown	25.9	NA	ABL1 Positive		
Well C	FAM	BCR-ABL1	Unknown	30.24	NA	BCR-ABL1 Positive	BCR-ABL1 Positive, ABL1 Positive	BCR-ABL1 Positive, ABL1 Positive
	ROX	ABL1	Unknown	26.34	NA	ABL1 Positive		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.