


Report Title

Test Name:	HPV-Q Comp	Test Sample:	4519052025
Test Time:	2025-05-19-16-34-14_1	Report Time:	2025-05-19

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	HPV18	Unknown	19.18	NA	HPV18 Positive	HPV18 Positive, HPV16 Positive, Other HPV Positive, IC Detected	HPV18 Positive, HPV16 Positive, Other HPV Positive, HPV6,11,8 Positive, HPV26,81, Positive, HPV40,43, Positive, IC Detected
	VIC	HPV16	Unknown	17.87	NA	HPV16 Positive		
	ROX	Other HPV	Unknown	19.06	NA	Other HPV Positive		
	CY5	IC	Unknown	18.23	NA	IC Detected		
Well B	FAM	HPV6,11,8	Unknown	19.0	NA	HPV6,11,8 Positive	HPV6,11,8 Positive, HPV26,81, Positive, HPV40,43, Positive, IC Detected	
	VIC	HPV26,81,	Unknown	20.28	NA	HPV26,81, Positive		
	ROX	HPV40,43,	Unknown	19.12	NA	HPV40,43, Positive		
	CY5	IC	Unknown	17.54	NA	IC Detected		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
 When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.