


# Report Title

Test Name:	HPV-Q	Test Sample:	HPVpc
Test Time:	2025-09-18-15-12-48_1	Report Time:	2025-09-18

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	HPV 18	Unknown	14.82	NA	HPV 18 Positive	HPV 18 Positive, HPV 16 Positive, HPV Other Positive, IC Detected	HPV 18 Positive, HPV 16 Positive, HPV Other Positive, IC Detected
	VIC	HPV 16	Unknown	14.56	NA	HPV 16 Positive		
	ROX	HPV Other	Unknown	13.91	NA	HPV Other Positive		
	CY5	IC	Unknown	12.25	NA	IC Detected		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;  
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.