


Report Title

Test Name:	Jak2-Q	Test Sample:	0611012025
Test Time:	2025-01-11-14-38-13_1	Report Time:	2025-01-11

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	WT	Unknown	20.3	NA	WT Positive	WT Positive, Mut Positive	WT Positive, Mut Positive
	VIC	Mut	Unknown	19.02	NA	Mut Positive		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.