


# Report Title

Test Name:	MTB-NTM	Test Sample:	222012202400
Test Time:	2024-12-20-15-37-06_1	Report Time:	2024-12-20

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	VIC	NTM	Unknown	21.46	NA	NTM Positive	NTM Positive, MTB Positive, IC Detected	NTM Positive, MTB Positive, IC Detected
	ROX	MTB	Unknown	19.67	NA	MTB Positive		
	CY5	IC	Unknown	18.28	NA	IC Detected		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;  
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.