


Report Title

Test Name:	MTC-Q	Test Sample:	MTCpc
Test Time:	2025-10-27-13-48-05_1	Report Time:	2025-10-27

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	IS6110	Unknown	23.15	NA	IS6110 Positive	IS6110 Positive, MPT64 Positive, IC Detected	IS6110 Positive, MPT64 Positive, IC Detected
	VIC	IC	Unknown	21.42	NA	IC Detected		
	ROX	MPT64	Unknown	22.39	NA	MPT64 Positive		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.