


# Report Title

Test Name:	PML-RARA	Test Sample:	4703062025
Test Time:	2025-06-03-13-24-57_1	Report Time:	2025-06-03

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	BCR-1	Unknown	13.45	NA	BCR-1 Positive	BCR-1 Positive, ABL1 Positive	BCR-1 Positive, ABL1 Positive, BCR-2 Positive, ABL1 Positive, BCR-3 Positive, ABL1 Positive
	ROX	ABL1	Unknown	14.76	NA	ABL1 Positive		
Well B	VIC	BCR-2	Unknown	14.31	NA	BCR-2 Positive	BCR-2 Positive, ABL1 Positive	
	ROX	ABL1	Unknown	14.75	NA	ABL1 Positive		
Well C	FAM	BCR-3	Unknown	15.56	NA	BCR-3 Positive	BCR-3 Positive, ABL1 Positive	
	ROX	ABL1	Unknown	14.84	NA	ABL1 Positive		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;  
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.