


# Report Title

Test Name:	SCA-Q	Test Sample:	1806062025
Test Time:	2025-06-06-13-31-19_1	Report Time:	2025-06-06

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	ROX	WT	Unknown	17.39	NA	WT Positive	WT Positive	WT Positive, Mut Positive
Well B	FAM	Mut	Unknown	19.68	NA	Mut Positive	Mut Positive	



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;  
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.