


Report Title

Test Name:	STI-Q Comp	Test Sample:	3006032025
Test Time:	2025-03-06-14-31-28_1	Report Time:	2025-03-06

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	CA	Unknown	20.97	NA	CA Positive	CA Positive, TG Positive, HD Positive	CA Positive, TG Positive, HD Positive, CA Positive, TG Positive, HD Positive, CA Positive, TG Positive, HD Positive
	VIC	TG	Unknown	20.61	NA	TG Positive		
	ROX	HD	Unknown	19.84	NA	HD Positive		
	CY5	IC	Unknown		NA			
Well B	FAM	CA	Unknown	24.03	NA	CA Positive	CA Positive, TG Positive, HD Positive	
	VIC	TG	Unknown	23.57	NA	TG Positive		
	ROX	HD	Unknown	23.01	NA	HD Positive		
	CY5	IC	Unknown		NA			
Well C	FAM	CA	Unknown	27.24	NA	CA Positive	CA Positive, TG Positive, HD Positive	
	VIC	TG	Unknown	26.71	NA	TG Positive		
	ROX	HD	Unknown	26.44	NA	HD Positive		
	CY5	IC	Unknown		NA			
Well D	FAM	CA	Unknown	32.82	NA	CA Positive	CA Positive, TG Positive, HD Positive	
	VIC	TG	Unknown	29.58	NA	TG Positive		
	ROX	HD	Unknown	29.64	NA	HD Positive		
	CY5	IC	Unknown		NA			



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.