



# Nucleic Acid Extraction Kits

---

- Extensive solution for various applications that ensure balance of high yield and reproducibility
- Delivers High Yield and Purity for Different Downstream Applications such as PCR, RT-PCR, Sequencing, NGS etc.

## TISSUE & BODY FLUID DNA Extraction Kit

- Optimized solution for cell/ tissue lyses & purification of intact genomic DNA from broad range of biological samples
- Can process sample using silica membrane based spin column/ magnetic beads technology

## Cell free DNA Extraction Kit

- Allows the isolation of cell-free DNA from either fresh or frozen plasma/serum samples
- Enables efficient purification of circulating nucleic acids from human plasma, serum, or urine

## FFPE DNA Extraction Kit

- Designed for the purification of genomic DNA from formalin- fixed paraffin-embedded (FFPE) tissue sections.
- Enables the partial reversing of the modification caused by formalin which results in high yield with good quality.

## SALIVA DNA Isolation Kit

- Isolation of genomic DNA from Human saliva samples and swabs
- Can process sample using silica membrane based spin column/ magnetic beads technology

## DBS DNA Extraction Kit

- Isolation of genomic DNA from dried blood spots
- Can process sample using silica membrane based spin column/ magnetic beads technology

## Fecal DNA Extraction Kit

- Isolate pure genomic DNA from fresh or frozen stool or feces samples
- Can process sample using silica membrane based spin column/ magnetic beads technology

## BLOOD DNA Extraction Kit

- Purification of genomic DNA from fresh whole blood/ blood stabilized with anticoagulants
- Unique buffer system for efficient RBC lyses and simultaneous recovery of WBCs by centrifugation

## Viral RNA Extraction Kit

- Available in both Pre-filled Ready-to Use Format and Bulk Packs
- Carrier RNA included to enhance quality of eluted viral RNA
- Recovery of  $\geq 95\%$  viral RNA

## Bacterial DNA/RNA Extraction Kit

- Allows isolation of gDNA/RNA from both gram negative and gram positive cultures
- Can process samples using silica membrane based spin column/ magnetic beads technology

## FORENSIC DNA Extraction Kit

- Isolate genomic DNA from cultured cells as well as various tissue types
- Can process sample using silica membrane based spin column/ magnetic beads technology

Commercial Name	Cat No.	Pack Size	Storage & Validity
SpiNXT Tissue & body fluids DNA Extraction Kit	G2M131420	50	RT, 12 Months
SpiNXT Tissue & body fluids DNA Extraction Kit	G2M131320	250	
MagNXT Tissue & body fluids DNA Extraction Kit	G2M212121	50	
MagNXT Tissue & body fluids DNA Extraction Kit	G2M212021	250	
MagNXT Tissue & body fluids DNA Extraction Kit	G2M211921	500	
MagNXT Tissue & body fluids DNA Extraction Kit	G2M211921-PF	480	
SpiNXT Saliva DNA Extraction Kit	G2M171721	50	
SpiNXT Saliva DNA Extraction Kit	G2M171821	250	
MagNXT Saliva DNA Extraction Kit	G2M191821	50	
MagNXT Saliva DNA Extraction Kit	G2M191921	250	
MagNXT Saliva DNA Extraction Kit	G2M192021-PF	480	
SpiNXT Forensic DNA Extraction Kit	G2M242421	50	
SpiNXT Forensic DNA Extraction Kit	G2M242521	250	
MagNXT Forensic DNA Extraction Kit	G2M242021	50	
MagNXT Forensic DNA Extraction Kit	G2M241921	250	
MagNXT Forensic DNA Extraction Kit	G2M241821-PF	480	
SpiNXT DBS DNA Extraction Kit	G2M232421	50	
SpiNXT DBS DNA Extraction Kit	G2M232321	250	
MagNXT DBS DNA Extraction Kit	G2M271921	50	
MagNXT DBS DNA Extraction Kit	G2M272021	250	
MagNXT DBS DNA Extraction Kit	G2M272121-PF	480	
SpinRNA Viral RNA Extraction Kit	G2M080820	250	
MagRNA-II Viral RNA Extraction Kit	G2M030420	500	
MagRNA-II Viral RNA Extraction Kit	G2M030520-PF	480	
MagRNA-II Viral RNA Extraction Kit	G2M030620	960	
SpiNXT Blood DNA Extraction Kit	G2M181820	50	
SpiNXT Blood DNA Extraction Kit	G2M181720	250	
MagNXT Blood DNA Extraction Kit	G2M181921	50	
MagNXT Blood DNA Extraction Kit	G2M182021	200	
MagNXT Blood DNA Extraction Kit	G2M182121	500	
MagNXT Blood DNA Extraction Kit	G2M182121-PF	480	
SpiNXT Bacterial DNA Extraction Kit	G2M141421	50	
SpiNXT Bacterial DNA Extraction Kit	G2M141521	250	
SpiNXT Bacterial RNA Extraction Kit	G2M121220	50	
SpiNXT Bacterial RNA Extraction Kit	G2M121320	250	
MagNXT Bacterial DNA Extraction Kit	G2M343421	50	
MagNXT Bacterial DNA Extraction Kit	G2M343521	250	
MagNXT Bacterial DNA Extraction Kit	G2M343621-PF	480	
MagNXT Bacterial RNA Extraction Kit	G2M282821	50	
MagNXT Bacterial RNA Extraction Kit	G2M282921	250	
MagNXT Bacterial RNA Extraction Kit	G2M283021-PF	480	
SpiNXT Cell free DNA Extraction Kit	G2M303021	50	
SpiNXT Cell free DNA Extraction Kit	G2M303121	250	
MagNXT Cell free DNA Extraction Kit	G2M313121	50	
MagNXT Cell free DNA Extraction Kit	G2M313221	250	
MagNXT Cell free DNA Extraction Kit	G2M313321-PF	480	
SpiNXT Fecal DNA Extraction Kit	G2M323421	50	
SpiNXT Fecal DNA Extraction Kit	G2M323521	250	
SpiNXT FFPE DNA Extraction Kit	G2M232121	50	
SpiNXT FFPE DNA Extraction Kit	G2M232221	250	



### Genes2Me Pvt. Ltd.

1105, 11th Floor, Tower B, SAS Tower, Medicity,  
Sector - 38, Gurgaon - 122 001, Haryana, India  
Tel : + 91 18001 214030 / + 91 88000 23600  
E-mail : contact@genes2me.com www.genes2me.com

